

INSTRUCTIONS FOR USE – ASK4MASK® MEDICAL MASK

I. Intended purpose and use

- The Type II R ASK4MASK® medical mask is intended to limit the transmission of infectious agents (viruses, bacteria, fungi) from healthcare personnel to patients during surgical procedures and in other medical environments with similar requirements.

II. Classification:

In accordance with SR EN 14683:2019+AC:2019:

- Type IIR: the ASK4MASK® medical mask has a bacterial filtration efficiency (BFE) $\geq 98\%$, breathability of $< 60 \text{ Pa/cm}^2$, microbial cleanliness $\leq 30 \text{ cfu/g}$, and splash resistance $\geq 16.0 \text{ kPa}$.

The ASK4MASK® medical mask is classified according to SR EN 14683:2019 +AC:2019, Chapters 4 and 5.2.7.

In accordance with Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009, and repealing Council Directives 90/385/EEC and 93/42/EEC, the ASK4MASK® medical mask is classified as a Class I medical device, in accordance with Annex VIII, Rule 1.

III. Instructions for use

These instructions apply both to medical masks packed in boxes of 5 pieces and those packed in boxes of 50 pieces.





EURO IMPACT S.R.L.

VAT No. RO 13881247 | Trade Register No.: J2001000283246

Registered office: 698 Ruscova Village, Maramureş County, PC 437260, Romania

Work point: 772 Leordina, Maramureş County, Romania

Phone: +40 783016964/+40 730336088, e-mail: office@ask4mask.ro,

web: <http://ask4mask.ro>

1. Before touching the mask, sanitise your hands using an alcohol-based disinfectant solution or wash them thoroughly with soap and water.
2. Check the mask for any tears or holes.
3. Identify the top of the mask. This is the side with the metal nose clip (rigid strip).
4. Identify the inside of the mask. This is the white side of the mask.
5. Ensure the coloured side is facing outward.
6. Place the mask over your face and position the ear loops behind your ears.
7. Pinch the metal nose clip (or rigid strip) to mould it to the shape of your nose.
8. Adjust the mask to cover your mouth and chin, ensuring there are no gaps between the face and the mask.
9. Do not touch the mask during use to avoid contamination.
10. If you accidentally touch the mask, sanitise your hands using an alcohol-based disinfectant or wash them thoroughly with soap and water.
11. To remove the mask, take off the ear loops without touching the front of the mask, keeping it away from your face. Keep the mask away from yourself and other surfaces to avoid contamination.
12. Dispose of the mask immediately in a closed bin and sanitise your hands.
13. Replace the mask as soon as it becomes damp.

In accordance with
the original

**CONFORM CU
ORIGINALUL**

Attention!

Do not reuse a used mask.

The ASK4MASK® medical mask is for single use only.



	Position	Name	Date	Signature
Prepared by	RMC	Marius Pop	02.12.2021	
Approved by	General Manager	Simion Pop	02.12.2021	